



# Immunofluorescence Has Never Been This Easy

#### ✓ Fast and Simple Handling

Simplify your staining procedure—perform all steps in one single slide


#### ✓ Cost-Effective Experiments


Reduce your costs—use only small numbers of cells and a low amount of medium and antibodies


### ✓ High-Resolution Imaging

Get brilliant microscopic images due to the slides' optical specifications









## **Saving Time With Immunofluorescence Assays**

### Comparison of Immunocytochemistry Protocols: Traditional Staining vs. Staining With ibidi Solutions


#### **Protocol With Cells on Coverslips**

Traditional method with nail polish mounting


- Sterilize coverslips and slides
- Coat coverslips
- Place sterile coverslips into 6-well plate
- Seed cells in large volume
- Peel off the coverslip
- Wash
- Fix wash permeabilize wash block
- Incubate in primary antibody wash incubate in secondary antibody
- \/\/ash
- Mount cells with mounting medium
- Mount coverslip with nail polish

#### Protocol With ibidi µ-Slides

Time-saving method using all-in-one chambers



- Sterilize coverslips and slides
- Coat coverslips
- Place sterile coverslips into 6-well plate
- Seed cells in large volume
- Peel off the coverslip
- Wash
- Fix wash permeabilize wash block
- Incubate in primary antibody wash incubate in secondary antibody
- Wash
- Mount cells with mounting medium
- Mount coverslip with nail polish



Fluorescence microscopy of rat fibroblasts (Rat1) in a µ-Slide 18 Well Glass Bottom. Red: alpha-tubulin; green: F-actin, stained with LifeAct-TagGFP2 Protein; blue: nuclei (ibidi Mounting Medium with DAPI). 60x objective lens, oil immersion.

#### ibidi Solutions for Immunofluorescence

#### Tailored for Your Assay: Choose from 3 Unique Solutions



#### **Chambered Coverslips**

- Up to 18 non-removable wells on a coverslip bottom
- Versatile use for different cell culture applications
- Different coatings available




#### **Channel Slides**


- Six parallel channels on a coverslip bottom
- Homogeneous cell and antibody distribution and small medium amounts
- Different channel heights and coatings available




#### **Chamber Slides**

- Removable silicone chambers on a standard glass slide
- Ideal for long-term storage
- Suitable for high-throughput screening







#### Immunofluorescence Has Never Been This Easy

ibidi provides several solutions that fit your needs for immunofluorescence assays:







Find many more varieties	Chambered Coverslips	Channel Slides	Chamber Slides, removable
on the ibidi website.	μ-Slide 2   4   8   18 Well	$\mu\text{-Slide VI}^{0.4} \mu\text{-Slide I Luer}$	3   8   12 Well Chamber
Bottom material	Glass Coverslip or Polymer Coverslip	Glass Coverslip or Polymer Coverslip	Standard glass slide
Microscope type	Inverted	Inverted	Inverted & upright
Mounting medium	Non-hardening	Non-hardening	Hardening
Sample storage	Short-term	Short-term	Long-term



### Support Your IF Assay With the Ready-to-Use ibidi Mounting Medium

- Available with and without DAPI
- Non-hardening—facilitiates the mounting of samples in channel slides
- Very low autofluorescence and prevention of photobleaching
- Allows sample storage for weeks without additional coverslips





#### Need a Detailed Guide?

Find more detailed information in our Application Guide.

#### Get Your Free Samples at: ibidi.com/free-samples

#### Front Page Image Information:

Top: Laser scanning microscopy of RDRGN and Schwann cells in a μ-Slide 8 Well, stained for neurofilament (green), NGFR (magenta), and DAPI (white). T. Weiss, Division of Plastic and Reconstructive Surgery, Medical University of Vienna, Austria.

Left: Trabecular meshwork cells of the human eye in the ibidi 8 Well Chamber, removable, stained for F-actin filaments (green) and DAPI (blue). Samantha Shan, School of Optometry, The Hong Kong Polytechnic University.

Middle: Epi-fluorescence of pMBMECs in the 12 Well Chamber, removable, stained for endothelial cell junctions (ZO-1, red) and DAPI (blue). S. Aydin, B. Engelhardt, R. Lyck, Theodor Kocher Institute, Bern, Switzerland.

Right: Widefield fluorescence of MDCK cells in the μ-Slide VI <sup>0.4</sup>, stained for F-actin (green), mitochondria (red), and DAPI (blue). Data by ibidi R&D.

ibidi GmbH Lochhamer Schlag 11 82166 Gräfelfing Germany Tel.: +49 89 / 520 46 17 - 0 Fax: +49 89 / 520 46 17 - 59

E-Mail: info@ibidi.de © ibidi GmbH, FL\_104, V 4.0, 2020/01 For more information on our products, please visit us at:

